




SEE CHART FOR DIMENSIONS



AMP VERTICALLY ADJUSTABLE HARDWARE

- Toggle lever allows the Amp pad to be positioned at multiple angles.



PRODUCT FEATURES

- + Amp support may be used with any wheelchair.
- + Installs onto either left or right side of seat cushion.
- + Toggle lever allows amp to be positioned at multiple angles.
- + Amputee pad has adduction to help center residual limb.
- + New foam construction has soft entry layer for comfort and firmer base layer for support.
- + Amputee pad is covered with incontinent proof fabric.

QUICK RELEASE AMP HARDWARE

- Fast and simple removal with the quick release assembly.



Attach to Cushion Bottom

10" AMP EXTENSION FLAT BOTTOM

model number	width x length
FAMP810	8 x 10
FAMP910	9 x 10
FAMP1010	10 x 10
FAMP1110	11 x 10
FAMP1210	12 x 10

10" AMP EXTENSION SLING BOTTOM

model number	width x length
RAMP810	8 x 10
RAMP910	9 x 10
RAMP1010	10 x 10
RAMP1110	11 x 10
RAMP1210	12 x 10

14" AMP EXTENSION FLAT BOTTOM

model number	width x length
FAMP814	8 x 14
FAMP914	9 x 14
FAMP1014	10 x 14
FAMP1114	11 x 14
FAMP1214	12 x 14

14" AMP EXTENSION SLING BOTTOM

model number	width x length
RAMP814	8 x 14
RAMP914	9 x 14
RAMP1014	10 x 14
RAMP1114	11 x 14
RAMP1214	12 x 14

NOTES

+ RIGID INSERT FOR SEAT CUSHION INCLUDED.

STANDARD COVER



COMFORT-TEK FABRIC: the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences.

Attach to Hard Seat Pan or for Bi-Lateral addition to FAMP/RAMP

10" EXTENSION KIT

model number	width x length
AMPKIT810	8 x 10
AMPKIT910	9 x 10
AMPKIT1010	10 x 10
AMPKIT1110	11 x 10
AMPKIT1210	12 x 10

14" EXTENSION KIT

model number	width x length
AMPKIT814	8 x 14
AMPKIT914	9 x 14
AMPKIT1014	10 x 14
AMPKIT1114	11 x 14
AMPKIT1214	12 x 14

NOTES

+ DOES NOT INCLUDE RIGID INSERT FOR SEAT CUSHION.