



HALO SAFETY DEVICE

HALO

SAFETY RING

The **ABC's** of Entrapment

- **Zone 1:** Within the rail.
- **Zone 2:** Between the top of the compressed mattress and the bottom of the rail, between rail ports.
- **Zone 3:** Between the rail and the mattress.
- **Zone 4:** Between the top of the compressed mattress and the bottom of the rail, at the end of the rail.
- **Zone 5:** Between the split bed rails.
- **Zone 6:** Between the end of the rail and the side edge of the head or foot board.
- **Zone 7:** Between the head or footboard and the mattress end.

*Mounts to Bed Frame*

Both Halo ALC and Institutional Bed Model mount to bed frame for stability.



FEATURES

The Halo Safety Ring, the only patented institutional bed mobility device that prevents lateral mattress movement, and gap space in the most critical zones: 1, 2, 3, & 4. The Halo Safety Ring meets and exceeds the HBSW and FDA's final dimensional guidelines against bed entrapment. The Halo Safety Ring passes the new "Bed System Measurement Device" pass/fail test that has been approved for Bed System testing by the FDA and state surveying teams.

- Provides exit points for neck, head, & chest and eliminates the four most critical zones of entrapment.
- Meets FDA guidelines for patient transfers & therapy bedside assistant for both patient transfers & therapy.
- Promotes bed mobility & increases environmental independence.
- Prevents lateral mattress movement and gap space elimination in the most critical zones: 1, 2, 3, and 4

BED HANDLE ROTATES

ROTATES 90° EITHER DIRECTION TO ACCOMMODATE DIFFERENT USER POSITIONS AND ALLOW EASE OF BED EXIT.



HALO SAFETY RING

INSTITUTIONAL BED MODEL



- Available in single and double sided for institutional style beds.
- Easy to install turn-buckle mounting bracket.
- 1000 lb. weight capacity steel construction

MODEL NUMBER	DESCRIPTION
77101	One sided circle device with mattress stay bar and hardware.
77102	Both sides circle device with hardware.
77103	One sided circle device with hardware for carrol.
77113	Halo Safety Ring, adjustable length bed (ALB) unit complete

HALO ALC

ASSISTED LIVING/HOME CARE MODEL



- The universal width and height adjustment accommodate all bed styles (up to queen) and allows bed to be off the wall, preventing bed roll-outs & falls.
- Comes standard with a universal height adjustment system to accommodate all bed systems' mattress and box spring heights.
- 500 lb. weight capacity steel construction.

MODEL NUMBER	DESCRIPTION
77121	Halo ALC (assisted living centers) complete kit

Halo Accessories

GAP SPACE ELIMINATOR


MODEL NUMBER - 77117



- Eliminates Zone 7 gap space.
- Designed to eliminate gaps between mattresses in cinrements of 2", 4", and 6".
- Units may be used individually or together, depending on size of gap to be eliminated.
- Each unit has a removable machine washable cover and a three strap system.

LEG PROTECTOR

MODEL NUMBER - 77116



- Transparent design eliminates any isolation for Dementia & Alzheimer's Patients.
- Prevents the upper from penetrating any one part of the Halo Circle.
- May be hand washed and air dried.

CIRCLE PAD PROTECTOR

MODEL NUMBER - 77118



- Prevents bruising from any medical disorder that might cause the skin to become thin, frail, or to bruise.
- Promotes increased safety and a better sleeping environment.
- Installs with simple hook and loop attachments.